

SPECIAL CONFIGURATIONS FOR THE CALCAP PROGRAM

Customizing the CALCAP Program

The CALCAP program can be customized in a variety of ways to accommodate different research protocols. The following features are available:

- C Summary output can be sent to the screen or to a printer. Printed output includes additional information explaining special codes and describing the normative comparison group.
- C At the end of each task the computer can give feedback on the subject's performance relative to age- and education-matched controls. Many subjects find this feedback helpful and motivating.
- C Entertaining music is played at the end of all choice reaction time measures. Some subjects enjoy this feature, others are annoyed by it. Music can be turned on or off at any time during the tasks.
- C The background can be either black or blue. All normative data is based on a blue background.
- C A demo version of the program can be invoked for instructional purposes or presentations.

Most of these features can be controlled from the command line when the CALCAP program is started. The different command line switches are described below and in the Figure on the next page.

Command Line Switches

Feedback. [/feed] The CALCAP program can provide feedback on performance at the end of each task. This feedback takes the form of statements designed to encourage or motivate the examinee ('That was very good,' 'You had a little trouble with that task—the next task will be quite different.') Each task has three levels of feedback: (1) performance above

expected levels; (2) performance in the average range; (3) performance below expected levels. Feedback is tailored to the age and education of the examinee. Default value is Feedback OFF. The /feed command line switch turns Feedback ON.

Music. [/nm] The CALCAP program plays a random 5-15 second musical selection at the end of each choice reaction time measure. Default value is Music ON. The /nm command line switch turns Music OFF.

Printer. [/p] At the end of the CALCAP test battery you are shown the test results on the video display. At this point you can also send the results to the printer. Alternatively, you can specify that you always want the results sent to the printer by using the /p option. Default value is Printer OFF. The /p command line switch turns Printer ON.

Text Resolution. [/novga] The original CALCAP program was designed for CGA monitors and all formal stimulus materials are displayed at CGA resolution. On EGA, VGA and XGA monitors, however, the CALCAP program will display text instructions using a high-resolution display font. You can force the CalcAP program to use CGA resolution for both instructions and stimulus materials by adding the /novga switch to the command line. Default is VGA Resolution ON. The /novga command line switch turns VGA Resolution OFF.

Multitasking Warning. [/nww] The CALCAP program needs exclusive use of the microprocessor in your computer. If it detects active multi-tasking software such as Microsoft® Windows®, it issues a warning. The /nww switch (No Windows Warning) eliminates this warning.

